

Press release

New Product Launch September 2009

Protect your staff and your instruments – introducing the versatile new clip trays from Prestige Medical

Prestige Medical say they have developed their new half size dental clip trays in order to give dentist the benefit of a seamless, hands-free process of instruments from the chair, through the cleaning and sterilization process - eliminating the risk of sharps injuries and damage.

Simple to use, the new trays are designed to safely carry the contaminated instruments from point of use to be loaded into the UltraClean II washer disinfector. Once the automatic cleaning process has been completed, the trays are removed and then placed directly into the C3 Advance autoclave. There is no need for transfer from one carrier to another, eliminating contact with the instruments and therefore the potential for sharps injuries.



Versatile new clip trays from Prestige

The new trays have been designed to take maximum advantage of the high capacity of the C3 Advance autoclave which enables it to carry a full load of 12 of the half size trays directly from the washer disinfector. This greatly speeds up the decontamination process and subsequently increases throughput.

Continued over ../

 **Prestige Medical**

Integrated decontamination solutions

Press release

The new clip trays are made in high quality grade 304 stainless steel and feature silicone moulded rests and retaining clips. They are supplied in sets of 2 trays (part number 309081) or 12 trays (part number 309082).

More information is available from Prestige Medical direct by calling 01254 844 103 or email to sales@prestigemedical.co.uk

For more information contact: Christine Bowness, Marketing Manager 01254 844 101
or email christine.bowness@prestigemedical.co.uk



Prestige Medical

Integrated decontamination solutions