

Press release

Prestige Medical launches new 'clean machine'

As part of their drive to provide dentists with a range of integrated decontamination solutions, Prestige Medical has announced the launch of UltraClean II —a bench top washer/disinfector/dryer which they claim will give guaranteed compliant performance to the 'best practice' standard in the HTM 01-05 guidelines.



(Caption: New UltraClean II washer/disinfector/dryer)

With a choice of 3 programme cycles —small load, normal load and a double wash for heavily contaminated instruments —UltraClean II is also claimed to have a high load capacity of 6 dental trays per cycle, making cleaning more cost effective.

When used in conjunction with an autoclave that can run in vacuum or non vacuum mode, Prestige Medical say that UltraClean II will allow a seamless, hands free process from the dentist, through the cleaning to sterilization in the autoclave, which should help to eliminate sharps injuries and instrument damage.

UltraClean II can be supplied with optional data logger or printer to enable data capture (necessary for compliance) and has full warranty and service support using their own team of service engineers.

More information is available from Prestige Medical direct by calling 01254 844 103 or email to sales@prestigemical.co.uk

For more information contact : Christine Bowness, Marketing Manager 01254 682 622
or email christine.bowness@prestigemical.co.uk



Integrated decontamination solutions